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# Problems with calibration of ear probes for measuring distortion product otoacoustic emissions

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**Abstract** *The intention of this article is to give a survey on ear probe calibration methods (Siegel 2002) and their fundamental problems. Moreover, the influence on the measurement of distortion product otoacoustic emissions (DPOAE) is described and ideas to improve calibration are presented.*

*The main goal of ear probe calibration is to generate a defined sound pressure level at the ear drum. Voltage at the loudspeaker is not functionally linked to sound pressure level at the eardrum due to the dependency of loudspeaker output level on the attached load impedance which differs individually. Main influencing factors are ear canal volume and ear drum impedance. Ear probe calibration is necessary to compensate for these impacts.*

*There are two major kinds of calibration methods (Whitehead et al. 1995). The most common calibration procedure is the in-the-ear adjustment strategy. Given a constant voltage at the loudspeaker, you measure the sound pressure at the ear probe microphone which is situated somewhere in the outer ear canal. The sound pressure level at the ear drum differs, dependent on the frequency, substantially from the measured sound pressure level at the microphone because of standing waves which result from the superposition of in-going and reflected waves. Especially at frequencies in the range of 3 to 7 kHz there*

could be fundamental deviations from the expected ear drum sound pressure level due to cancellations at the microphone place. These differences could exceed 20 dB SPL. Therefore, network models are useful to compensate for these effects. These models predict the differences between the sound pressure level at the microphone and at the ear drum on the basis of measurements of certain parameters as ear canal length and ear drum impedance. Indeed, the estimation of difference is very difficult, because the measurement of the required parameters is due to several side effects in the ear quite complex and thus imprecise. All in all, no complete correction will be possible.

Another calibration method is the constant voltage calibration, which measures the transfer function within a reference coupler and not within the ear. The sound pressure level is measured at the coupler microphone given a defined input voltage at the ear probe loudspeaker. This curve is used as reference transfer function and is shifted according to the individual ear canal volume. The advantage over the in-the-ear calibration is that the coupler microphone is located at the ear drum position, that means at the end of the coupler, so that there are no effects due to standing wave minima. Sources of error are deviations between the individual ear canal and the coupler transfer function. Further on,  $\lambda/2$ -resonances which are measured in the ear coupler could deviate from resonance frequencies in a real ear. Also, the amplitude of the resonance could differ enormously between coupler and ear.

The main advantage of this calibration strategy is the irrelevance of standing wave minima at the ear probe microphone place, so that the difference between assumed and real sound pressure level at the ear drum should be lower than with in-the-ear calibration in a frequency range up to about 5 kHz.

Considering all the additional undesirable effects in real ears, which were described above and which falsify the sound pressure level at the ear drum depending on the calibration method, it is apparent that calibration has a strong influence on DPOAE measurements, which are dependent on the primary tone level setting. The two primary tone levels  $L_1$  and  $L_2$  could deviate differently from their demand value, because they are located at different frequencies. Thus the optimal DPOAE amplitude, resulting from an ideal stimulus level setting (scissor paradigm:  $L_1 = 0,4L_2 + 39$  – Janssen et al. 1995; Kummer et al. 2000) will be changed. This results in a calibration-dependent varying course of the DPOAE growth function. Therefore, the calibration strategy has to be improved in order to assure a reliable acquisition and evaluation of DPOAE data, which is important for further trustworthy diagnostic conclusions.

To improve the in-the-ear calibration strategy it is necessary to develop reliable compensation models and test their quality in real ears. Further on, constant voltage calibration could possibly be enhanced by using various reference curves and chose one of these dependent on ear canal length measurement. Moreover, it is important to think of methods for evaluating the quality of a particular calibration strategy when used in real ears, since it is not easily possible to measure the sound pressure level at the ear drum. Subjective comparative tests could be helpful to solve this problem.

**Key words**

calibration  
in-the-ear calibration  
constant voltage calibration  
standing wave  
eardrum impedance  
ear canal length  
distortion product otoacoustic emissions

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
Fig. 1: Schematical presentation of the in-the-ear calibration procedure – left: measurement setup for the in-the-ear calibration – ear probe (2 loudspeakers, 1 microphone) placed in the outer ear canal; right, at the top: time signal presented at the loudspeaker (frequency modulated sine signal) – right, at the bottom: amplified time signal recorded at the microphone (superposition of the input signal and the signal reflected from the eardrum)

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Fig. 2: Sound pressure level at the ear probe microphone related to the loudspeaker input voltage using an ER-10C ear probe. The measurement was achieved in a syringe with a distance from the ending of the syringe to the ear probe microphone of 2 cm.

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Fig. 3: Illustration of the synthesis of standing waves – dark grey line: sine signal transmitted by the acoustic source; light grey line: completely reflected sine signal at a rigid ending; black line: resulting superposition of the two sound waves (standing wave). The points mark the locations where there are nodes (constant zero points) of the standing wave. The maxima of the standing wave are located at  $0$ ,  $\lambda/2$  and  $\lambda$  for the plotted range.

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Fig. 4: Absolute value and phase of the eardrum impedance according to a model by Shaw and Stinson (Shaw and Stinson 1981, compare Letens 1988).

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Fig. 5: Standing wave at a partly reflecting ending ( $Z < \infty \rightarrow$  amplitude of the reflected wave: 80%, phase shift:  $90^\circ$ ) – dark grey line: sine signal transmitted by the acoustic source; light grey line: partly reflected sine signal at a non-rigid ending; black line: resulting superposition of the two sound waves (standing wave). The vertical line marks the first maximum at a distance to the ending smaller than  $\lambda/2$ . The points mark the now not any more existing nodes at  $\lambda/4$  and  $3\lambda/4$ .

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
Fig. 6: User interface for a partly automatic ascertainment of the individual ear canal length by detecting the minimum of the standing wave. Adjustable parameters are the assumed model of the terminating impedance of the ear canal (rigid ending, model by Shaw/Stinson, and by Hudde/Engel), the number of averages for the measurement, the temperature, and the window width for the local minimum detection window (Fröstl 2001).

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Fig. 7: Equivalent circuit diagram for the acoustical impedance measurement –  $q_0$ : sound energy flux,  $Y_0$ : source admittance,  $p$ : sound pressure measured at the ear probe microphone,  $Y_x$ : load admittance (eardrum, ear canal) (Lodwig 1997).

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Fig. 8: Sound pressure level at the ear simulator (Brüel & Kjær 4157) microphone at a given voltage of 0,1 V at the ear probe (ER-10C) loudspeaker.

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Fig. 9: Schematic comparison of the deviation between actual and target sound pressure level for in-the-ear and constant-voltage calibration (Whitehead 1995).

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Fig. 10: Measurement of DPOAE growth functions at 4 kHz with constant-voltage calibration (black) and in-the-ear calibration (grey) at a subject with an estimated ear canal length of 22 mm ( $\rightarrow$  Minimum of the standing wave is located around 4000 Hz) – left: original DPOAE growth functions; right: DPOAE growth functions with the function measured with in-the-ear-calibration shifted 13 dB SPL along the  $L_2$  axis.

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- Janssen T, Kummer P, Arnold W (1995) Wachstumsverhalten der Distorsionsproduktemissionen bei kochleäre Hörstörungen. *Otorhinolaryngol Nova* 5: 34–46
- Keefe DH, Ling R, Bulen JC (1992) Method to measure acoustic impedance and reflection coefficient. *J Acoust Soc Am* 91 (1): 470–485
- Kummer P, Janssen T, Hulin P, Arnold W (2000) Optimal  $L_1$ - $L_2$  primary tone level separation remains independent of test frequency in humans. *Hearing Res* 146 (1): 47–56
- Letens U (1988) Über die Interpretation von Impedanzmessungen im Gehörgang anhand von Mittelohr-Modellen. Dissertation, Ruhr-Universität Bochum, Fakultät für Elektrotechnik
- Lodwig A (1997) Ein System zur Anpassung von Hörgeräten: Messtechnik, patientenbezogene Modellierung, Optimierung der Komponenten. Dissertation, Ruhr-Universität Bochum, Fakultät für Elektrotechnik
- Neely ST, Gorga MP (1998) Comparison between intensity and pressure as measures of sound level in the ear canal. *J Acoust Soc Am* 104 (5): 2925–2934.
- Shaw EAG, Stinson MR (1981) Network concepts and energy flow in the human middle-ear. *J Acoust Soc Am, Supplement* 1, 69: 43
- Siegel JH (1994) Ear-canal standing waves and high-frequency sound calibration using otoacoustic emission probes. *J Acoust Soc Am* 95 (5): 2589–2597
- Siegel JH, Dreisbach LE (1995) Optical placement of a probe tube in the occluded human ear canal. *Abstr Assoc Res Otolaryngol* 18: 119
- Siegel JH, Hirohata ET (1994) Sound calibration and distortion product otoacoustic emissions at high frequencies. *Hearing Res* 80 (2): 146–152
- Siegel JH (2002) Calibrating Otoacoustic Emission Probes. In: Robinette MS, Glattke TJ (eds.) *Otoacoustic Emissions – Clinical Applications* (second edition). Thieme Verlag, New York: 416–441
- Whitehead ML, Stagner BB, Lonsbury-Martin BL, Martin GK (1995) Effects of ear-canal standing waves on measurements of distortion-product otoacoustic emissions. *J Acoust Soc Am* 98 (6): 3200–3214

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*Dipl.-Ing. Jörg Müller: Born 1976. Studies in electrical engineering and information technologies between 1997 and 2002 at the Technische Universität München. Since 2002 academic employee at the experimental audiology laboratory of the ENT department at the Klinikum Rechts der Isar in Munich. Main research activity: distortion product otoacoustic emissions and their applicability in audiological diagnostics.*